


**Amptek Inc.**  
**14 Deangelo Drive**  
**Bedford, MA 01730 USA**  
**[www.amptek.com](http://www.amptek.com)**

# **Overview of XRF Instrumentation:**


## **Software**

### **AMPTEK, INC.**

14 Deangelo Drive, Bedford, MA 01730  
Ph: +1 781 275 2242 Fax: +1 781 275 3470  
[sales@amptek.com](mailto:sales@amptek.com) [www.amptek.com](http://www.amptek.com)



- Excitation source
  - X-ray tube or radioisotope
- Spectrometer
  - X-ray detector
  - Signal processing electronics
- Software
  - Spectrum correction and processing software
- Other
  - Radiation shielding
  - Sample fixture



# Software

AmpTek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)


## What does the software do?

### – Acquisition software

- Configures the hardware
- Starts and stops data acquisition
- Reads the raw pulse height spectrum from spectrometer

### – Analysis software



- Applies energy calibration
- Corrects raw spectrum for various physical effects and artifacts
- Identifies photopeaks
- Computes intensity (counts/sec) of each photopeak
  - Subtracts background
  - Subtracts counts due to overlapping peaks
- Computes concentrations
  - Corrects for matrix effects
  - Applies calibration factors



# Software


Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## What the software does, part 1: Acquires the raw spectrum



Raw spectra of stainless steel: linear (left) shows the main peaks while log plot (right) shows background and artifacts


Measured using Amptek Experimenter's Kit, Ag anode tube, W/Al filter, X-123-SDD spectrometer, and XRS-FP analysis software




# Software

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## What the software does, part 2: Processes the spectrum



After background removal




After peak deconvolution

Corrects for physical effects (pileup, escape peaks, scatter)

Removes background continuum

Fits to sum of peaks



# Software

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

## What the software does, part 3: Quantifies results

Element	Intensity (cts/sec)		
V	0.8	±	0.8
Cr	359.8	±	7.3
Mn	16.2	±	2.1
Fe	975.2	±	12.0
Co	0.0	±	1.2
Ni	128.6	±	4.4
Cu	2.6	±	0.9
Mo	461.1	±	8.3

Peak intensities


Element	Concentration		
V	0.06%	±	0.07%
Cr	19.20%	±	0.39%
Mn	0.81%	±	0.10%
Fe	64.66%	±	0.80%
Co	0.00%	±	0.00%
Ni	12.05%	±	0.42%
Cu	0.17%	±	0.06%
Mo	3.05%	±	0.05%

Elemental concentrations

Deconvolution → intensity (X-rays/sec) in each peak

Applies calibrations and matrix corrections

Computes elemental concentrations (final answer)




# Software

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

## Types of XRF analyses


- Qualitative analysis
  - Identifies which elements are present (above a detection limit)
- Semi-quantitative analysis
  - Measures elements concentration using a Fundamental Parameters (FP) method to correct for matrix effects
  - *Standardless FP*: Algorithm estimates calibration and matrix factors
  - *Hybrid FP*: Use mix of calibration and algorithm
  - Relative accuracy of 5-10% (hybrid) to 10-20% (standardless)
- Quantitative analysis
  - Strict reference method: X-ray intensity is compared to a calibration curve from standard reference materials
  - Calibration must use same matrix
- All methods apply equally to bulk concentration, layer density ( $\text{mg}/\text{cm}^2$ ), layer thickness, etc.



**Amptek Inc.**  
**14 Deangelo Drive**  
**Bedford, MA 01730 USA**  
**[www.amptek.com](http://www.amptek.com)**

# Energy Calibration







# Energy Calibration

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## What is energy calibration?

- The spectrometer outputs a "pulse height spectrum", a histogram of the number of X-rays in each pulse height bin or "channel"
- Spectrometer does not "know" the energy associated with each pulse height bin; this can depend on many instrument parameters
- Energy calibration finds relationship between pulse height and energy: the "gain" (eV/channel) and "offset" (generally not zero)







# Energy Calibration

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## What is energy calibration?

- The spectrometer outputs a "pulse height spectrum", a histogram of the number of X-rays in each pulse height bin or "channel"
- Spectrometer does not "know" the energy associated with each pulse height bin; this can depend on many instrument parameters
- Energy calibration finds relationship between pulse height and energy: the "gain" (eV/channel) and "offset" (generally not zero)







# Energy Calibration

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## How is an energy calibration done?


- In a turn-key analyzer, usually done by manufacturer but user may need to check
- Start with a known spectrum, and do linear fit between known energy and channel number
- Should be at least two peaks which are (a) well above background, (b) separated from overlapping peaks, and (c) far apart in energy





**Amptek Inc.**  
**14 Deangelo Drive**  
**Bedford, MA 01730 USA**  
**[www.amptek.com](http://www.amptek.com)**

# Photopeak spectra






# Photopeaks

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## What are the energies of the photopeaks?

- Moseley's law gives approximate values; exact values are given in various tables
- Energy of each X-ray peak equals difference between atomic energy levels






# Photopeaks


Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)


## Example spectra

Tin spectrum




Multielement spectrum





**Amptek Inc.**  
**14 Deangelo Drive**  
**Bedford, MA 01730 USA**  
**[www.amptek.com](http://www.amptek.com)**

# Spectrum Processing




# Spectrum Processing

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

## What does the processing do?

- There are features in the spectrum, in addition to photopeaks from sample, which must be addressed for accurate analysis
  - Escape peaks
  - Escape continuum
  - Compton and Rayleigh scatter peaks
  - Brehmstrahlung continuum
  - Sum peaks
  - Overlapping peaks (a.k.a. interferences)
- Some of these cause photopeak counts to occur in other energy channels. The software adds these "lost" counts back into the photopeak
- Some of these can obscure real peaks. The software "cleans" the spectrum so real peaks are clearer







# Spectrum Processing

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## Escape peaks

- When an X-ray stops in a detector, by photoabsorption, a silicon atom in the detectors is in an excited state, so emits Si characteristic X-rays
- Some Si  $K_{\alpha}$  X-rays escape the detector  $\rightarrow$  false peak 1.75 keV below primary peak
- Must be corrected in software: reassigns lost counts to photopeak and away from background region







# Spectrum Processing

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## Background

- Several processes produce a background continuum:
  - Brehmstrahlung X-rays from tube scatter from sample or target
  - Photoelectrons escaping detector form a continuum
  - Continuum is mostly smooth but has steps, curves, etc – complex shapes
- Background obscures real peaks and degrades measurement precision/accuracy
- Software must distinguish background counts from real counts in same region







# Spectrum Processing

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## Sum peaks

- If two X-rays interact in detector at same time, detector sees one pulse with sum of energies
- Most important for the major peaks and at the highest count rates
- The sum peaks can obscure real peaks and counts there are lost from main peak







# Spectrum Processing

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## Overlapping peaks

- If another peak overlaps the spectral peak (region of interest) of the element to be analyzed, this is termed an "interference"
- Overlapping peaks can arise from sample itself (e.g. Fe  $K_{\alpha}$  peak overlaps Mn  $K_{\beta}$  peak) or from the environment (air contains Ar so measurement in air show Ar  $K_{\alpha}$ )
- Software must account for the overlap to obtain the correct intensity





# Spectrum Processing

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)


## How does one determine the area of the photopeak?

Several different methods can be used

- Direct sum of counts within a region of interest (ROI)
- Fit photopeak shape to a Gaussian, using a predetermined centroid and width
- Fit photopeak shape to a Gaussian, adjusting centroid and width during fit
- Fit photopeak shape to a more complicated theoretical curve (e.g. to include a tail or other features)
- Fit photopeak shape to a measured curve (reference spectrum)


Which is best?

- Simpler ones are quick and work well for large, separated peaks
- More complicated ones help with overlapping peaks, complex spectra, and intensities near the detection limits



**Amptek Inc.**  
**14 Deangelo Drive**  
**Bedford, MA 01730 USA**  
**[www.amptek.com](http://www.amptek.com)**

# Quantification





# Quantification

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com


## Calibration factors

- Intensity of each element's X-ray peak is proportional to the number of the atoms excited, i.e. concentration
- Only true if **everything** is held constant.
- Calibration factors depend on geometry, excitation energy, absorption in sample, many other factors
- An error in calibration → poor accuracy




## Matrix effects

- Every atom in the sample affects the intensity of X-rays from the analyte
- Attenuation
  - X-rays produced in the sample are attenuated by the sample
  - Attenuation depends on composition of sample
- Enhancement
  - If you add a high-Z element to a sample, it will excite the lower Z elements, thus increasing the intensity
- You must know the composition of the matrix to correct matrix effects







# Quantification

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)


## How to approach calibration & matrix effects

### – Strict "quantitative" analysis

- A reference based method, requiring standards with compositions that are well known and similar to those of material to be analyzed
- Standards are analyzed, intensities obtained, and a calibration plot is generated for elements of interest

### – Comparative analysis

- In many practical applications, one does not need concentrations per se but can compare to a reference
- In manufacturing, one often is often comparing items to predetermined reference samples which set limits



# Quantification

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)


## How to approach calibration & matrix effects

### – Standardless fundamental parameters (FP) analysis

- Theory is used to compute calibration factors and matrix corrections. Based on the "Sherman Equation"
- No calibration needed, so very easy. One does not need to obtain reference samples similar to what will be measured
- Equations are approximate, does not include real world tolerances, so relative accuracy is limited to ~10% best case

### – Hybrid FP or FP with standards

- Calibration is used to determine some factors while others, e.g. matrix corrections, are derived from theory



# Quantification

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

## Measurement equations

*Simplified but illustrative*

### – Equation [1]


- $C_i$  is the concentration of element  $i$
- $I_i$  is the X-ray intensity of element  $i$  (in counts/sec)
- $K_i$  is the calibration coefficient
- $M$  is a matrix correction factor
  - The “matrix” is all the other elements in the sample besides the analyte
  - X-rays emitted by every element in the sample interact with the atoms from all the other elements

$$C_i = I_i \cdot K_i \cdot M \cdot S \quad [1]$$

$$I_i = (G_i - B_i - P_i + \gamma_i) / T \quad [2]$$


### – Equation [2]

- $G_i$  is the gross or total counts in a region of interest (ROI) for the peak of element  $i$
- $B_i$  is the background counts
- $P_i$  is counts from overlapping or interfering peaks
- $\gamma_i$  is correction for artifacts
- $T$  is the measurement time



## Measurement uncertainty

- Every term in equations [1] and [2] contributes to measurement uncertainty
- There will be random errors (affecting precision)
  - Dominant source is usually from random counting statistics
  - If there are 100 counts in photopeak, relative precision will be 10%
  - To get 0.1% precision you need 1 million counts in peak
- There will be systematic errors if anything changes from calibration/model
  - If sample geometry changes, then calibration factor changes
  - If you change  $T_{\text{peak}}$ , HV on tube, filter on tube, etc then errors results
- There will be systematic errors due to approximations in models
  - Photopeaks are never truly Gaussian. A Gaussian deconvolution is only approximate
  - Escape peak models, sum peak models, background models only approximate
  - These errors are largest when a small peak of interest is on a large peak to be modeled



# Software


Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

## Measurement uncertainty

- Expanded uncertainty

$$\frac{\sigma_C^2}{C^2} = \frac{\sigma_K^2}{K^2} + \frac{\sigma_M^2}{M^2} + \frac{\sigma_S^2}{S^2} + \frac{(\sigma_G^2 + \sigma_B^2 + \sigma_P^2 + \sigma_\gamma^2)}{I^2}$$

- Each term includes random errors (due to counting statistics), systematic errors due to changing conditions (i.e. geometry), and systematic errors due to approximations in the models



# Amptek Software



Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)


## DPPMCA

- Data acquisition & control
- Limited analysis
  - Energy calibration
  - Define ROIs, compute net counts
- Free with processor

## Software developer kit

- User can integrate Amptek DPP with their analysis & display
- Cross-platform available for non-Windows





# Amptek Software

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

## XRS-FP

- Full analysis software
  - Processes spectra, fits photopeaks to get intensity, applies matrix corrections & calibration factors
- Many options
  - Standardless FP, FP with standards, purely standards based
  - Multiple excitation conditions, scattered intensity corrections, multiple thin film analyses
  - Gaussian peaks, ROI, reference spectra
- Can be customized for specific applications or user interfaces

