


**Amptek Inc.**  
**14 Deangelo Drive**  
**Bedford, MA 01730 USA**  
**[www.amptek.com](http://www.amptek.com)**

# **Overview of X-Ray Fluorescence Analysis**

## **AMPTEK, INC.**


14 Deangelo Drive, Bedford, MA 01730  
Ph: +1 781 275 2242 Fax: +1 781 275 3470  
[sales@amptek.com](mailto:sales@amptek.com) [www.amptek.com](http://www.amptek.com)




Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

# What is X-Ray Fluorescence (XRF)?

- A physical process:  
Emission of characteristic "secondary" (or fluorescent) X-rays from a material that has been excited by high-energy X-rays or gamma rays.
- A technique in analytical chemistry:  
Method to identify elements in a sample and measure their concentrations  
  
Non-destructive, quick, and simple to carry out.




- Incoming radiation hits an atom
- Ejects an electron from an inner shell, creating a vacancy
- An electron from an outer shell “drops down” to fill the vacancy.



- The excited atom emits an X-ray with energy equal to the difference between the levels


$$E_{Xray} = \Delta E = E_K - E_L$$

- Since each element has a unique set of levels, it produces a unique set of “characteristic” X-rays




# Physical Process

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com




		Fe	Ni	Cu	Zn	Pb
Atomic Levels						
K	1s	7,112	8,333	8,979	9,659	88,005
L <sub>1</sub>	2s	845	1,009	1,097	1,196	15,861
L <sub>2</sub>	2p <sub>1/2</sub>	720	870	952	1,045	15,200
L <sub>3</sub>	2p <sub>3/2</sub>	707	853	933	1,022	13,035
M <sub>1</sub>	3s	91	111	122	140	3,851
M <sub>2</sub>	3p <sub>1/2</sub>	53	68	77	91	3,554
M <sub>3</sub>	3p <sub>3/2</sub>	53	66	75	89	3,066
Characteristic X-Ray Lines						
K <sub>α1</sub>	K – L <sub>3</sub>	6,404	7,478	8,048	8,639	74,969
K <sub>α2</sub>	K – L <sub>2</sub>	6,391	7,461	8,028	8,616	72,804
K <sub>β1</sub>	K – M <sub>3</sub>	7,058	8,265	8,905	8,572	84,936
L <sub>α1</sub>	L <sub>3</sub> – M <sub>5</sub>	705	852	930	1,012	10,552
L <sub>α2</sub>	L <sub>3</sub> – M <sub>4</sub>	705	852	930	1,012	10,450
L <sub>β1</sub>	L <sub>2</sub> – M <sub>4</sub>	718	869	950	1,035	12,614


- XRF is similar to optical spectroscopy but at higher energy
- Independent of chemical state → **Elemental analysis**



- Intensity of X-ray line proportional to number of atoms → **Quantitative**
- X-rays pass through surface into sample
  - **Nondestructive** and **no sample preparation** is necessary
    - Best accuracy requires sample preparation
  - **Bulk measurement** rather than only surface
    - Notion of “bulk” vs “surface” depends on the X-ray energy




- Presence of Cu and Zn K lines → Elements are in sample (qualitative)
- Intensity of the lines → How much is in sample (quantitative)



## Typical spectrum and results

- Photopeak intensity varies with concentration
- Final result is quantitative concentration

		Certified		XRF	
High	Cr	1000 ± 20		895 ± 198	
	Br	1100 ± 22		1089 ± 23	
	Cd	300 ± 6		264 ± 28	
	Hg	1100 ± 22		1050 ± 53	
	Pb	1200 ± 24		1184 ± 39	
Low	Cr	401 ± 8		388 ± 167	
	Br	500 ± 10		487 ± 13	
	Cd	100 ± 5		68 ± 13	
	Hg	200 ± 5		183 ± 27	
	Pb	400 ± 8		398 ± 23	
Blank	Cr	0 ± 5		7 ± 40	
	Br	0 ± 5		1 ± 2	
	Cd	0 ± 5		9 ± 10	
	Hg	0 ± 5		0 ± 0	
	Pb	0 ± 5		10 ± 9	



# Measurement

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

## What are the main factors limiting XRF?


### – Counting variance

- The measurement is based on count discrete X-rays
- Arise from random processes → Inherent statistical variation in number of X-rays
- Percent uncertainty =  $1/\sqrt{N}$ 
  - 100 X-rays detected → 10% precision
  - 1,000 X-rays detected → 3% precision
  - 1,000,000 X-rays → 1000 ppm precision
- Good precision means many X-rays which means high count rates or long times

### – Detector response

- Photopeak has some width
- There is always spectral background and overlapping peaks
- Ability to remove these depends on counting variance, energy resolution, and accuracy of software algorithms
- Better energy resolution helps but there are physical and practical limits






# Measurement


Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

## What are the main factors limiting XRF?

### – Attenuation lengths

- Penetration depth depends on energy & therefore element
  - In silica, Al X-rays go 3  $\mu\text{m}$  while Sn go 3 mm
- Response depends on energy/element
- Sample condition & homogeneity are critical





# Measurement FAQ

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

- **How accurate is EDXRF?**


- In the best case, relative accuracy  $\sim 0.2\%$  ( $1.00\% \pm 20$  ppm).  
Requires sample prep, a known matrix, good statistics, etc
- Nondestructive screening, relative accuracy  $\sim 2\%$  ( $1.00 \pm 0.02\%$ )  
Requires careful optimize and setup, known sample type
- Quick check on unknown, relative accuracy  $\sim 20\%$  ( $1.0 \pm 0.2\%$ )

- **What is the detection limit for EDXRF?**

- $<1$  ppm for prepared samples in a known matrix under good conditions
- 10 ppm in nondestructive screening with no interfering elements
- When elements interference or overlap, 1% of other element

- **What elements can be analyzed with EDXRF?**

- Na to U (down to Be with EDS)
- Low Z elements (below S) are a challenge
- Need multiple measurements to cover a wide range of elements



# Measurement

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

## **XRF is one of many methods used in material analysis**

### **Advantages of XRF**


- Non-destructive
- No sample preparation
- Fast (seconds to minutes)
- Good precision and accuracy
- Measure Na to U
- Suitable for portable equipment and field use

### **Disadvantages of XRF**

- Limits of detection modest (10 ppm typical)
- Accuracy usually modest (few % relative)
- Difficult to use for lower Z elements

### **Best results require**


- Sample preparation (damaging)
- System optimization
- Matched calibration standards



Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

# XRF Applications




How is XRF used?




# Hazardous Material Screening

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

- Is there cadmium on this toy?
- Is there lead in this paint?
- Does this circuit board contain Pb, Cd, or Cr?
- **Nondestructive critical for screening products!**











# Metal Alloy Analysis

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com


- Are these bolts stainless steel 316?
- Is there Ni in this scrap metal?
- Is there Cd plating on this MILSPEC connector?
- How much Au is in a white gold ring?
- **Speed critical, accuracy moderate**



- How did the artist make their paints?
- Is this an ancient or a modern pigment?
- What is the effect of cleaning on the surface of a statue?
- **Nondestructive testing is vital for art!**










# Process Control

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com

- Is there any change in spectrum?
- Absolute composition not needed but quick, real-time, non-destructive vital.








# Field Measurements

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
www.amptek.com


Measuring ores in mines.




Identifying minerals on Mars.



Measuring soil contamination







**Amptek Inc.**  
**14 Deangelo Drive**  
**Bedford, MA 01730 USA**  
**[www.amptek.com](http://www.amptek.com)**


# Related Analytical Methods

## Energy Dispersive X-ray Spectroscopy

– a.k.a. EDS, EDX, XEDS, EDXA

- Uses electron beam in vacuum chamber to excite the atoms
- Electrons have short range in matter → Only way to measure lightest elements, down to Be (Z of 4)
- High spatial resolution (75  $\mu\text{m}$  spatial, 1  $\mu\text{m}$  in depth)





## Related

Amptek Inc.  
14 Deangelo Drive  
Bedford, MA 01730 USA  
[www.amptek.com](http://www.amptek.com)

# Wavelength Dispersive X-ray Spectroscopy (WDXRF)

- Crystal diffractometer disperses the X-ray wavelengths much like a prism disperses visible light.
- X-rays at a particular wavelength (energy) are recorded by a detector.
- It measures only one energy at a time; It obtains a spectrum by sweeping the wavelength over time
- Advantages of WDXRF
  - Much better energy resolution
  - Leads to much better accuracy and detection limits
- Disadvantage of WDXRF
  - Very long time to acquire whole spectrum
  - Requires destructive sample preparation
- Uses similar detectors and signal processors